

ABSORBENT COTTON WOOL (BPC)

DISCRIPTION:

It consists of good trichomes or good quality new comber obtained from seed coats .it is white and particularly odorless and composed of carefully carded fibers of average staple length.it offers appreciable resistance when pulled and does not shed a significant quantity of dust shaken gently.

FEATURES:

- High fluid and blood absorbency power.
- Chemically inert and soft to give maximum protection and not cause irritation.
- Free of compensatory coloring matters, foreign fibers and napes and external impurities.

INDICATIONS:

- Used for cleaning wounds
- For dampening and using for cleaning skin and wounds
- For padding and protection for particularly sensitive skin
- Use in surgical procedures for the efficient absorption of fluid and blood.

PRECUATIONS:

- Not really for applying onto open wounds as the free fibers can break away and stay in the wound.
- Sterile before use.

STORAGE: keep in cool & dry place.

Its absorbency may be reduce considerably by medications, prolong storage or by exposure to heat.

AVAILABLE SIZES:

- Absorbent cotton wool50gm
- Absorbent cotton wool100gm
- Absorbent cotton wool200gm
- Absorbent cotton wool250gm
- Absorbent cotton wool300gm
- Absorbent cotton wool400gm
- Absorbent cotton wool500gm
- Cotton ball 100gm
- Zig zag 100gm
- Orthoban cotton 5"



COTTON CRAFT (PVT) LTD PAKISTAN

Plot #407, 408- Sundar Industrial Estate, Raiwind Road, Lahore
Ph : +92-42-35298604-06 , Email: info@cottoncraftpvt.com

Website : www.cottoncraftpvt.com