

ABSORBENT GAUZE (BPC)

DISCRIPTION:

Cotton cloth of plain weave bleached to a good white, clean, and reasonably free from weaving defects, cotton leaf and shells.it is colorless or almost odorless.it may be slightly off-white if it has been sterilized.

FEATURES:

- Highly absorbent 100% cotton woven fabric.
- High fluid and blood absorbency power.
- Free of compensatory, foreign fibers and external impurities and weaving defects.
- Chemically inert and soft to give maximum protection and not cause irritation to skin.

INDICATIONS:

- Used to stop the flow of blood absorb drainage.
- Cushion the injured area.
- Provide a safeguard against contamination.
- Hold a medicated dressing in place and hold a splint in position.
- Immobilize an injured part of the body to prevent further injury and to facilitate healing.

PRECUATIONS:

- Sterile before use.
- Do not use if sterile packet is torn or ruptured.

STORAGE:

keep in cool & dry place.
Its absorbency may reduce considerably by medications, prolong storage.

AVAILABLE SIZES:

- Absorbent Gauze 100cmx100cm
- Absorbent Gauze 1x2.5M (ST+non-ST)
- Absorbent Gauze 1x10 M
- Absorbent Gauze 1x20M
- Absorbent Gauze 1x30M
- Absorbent Gauze 1x40M



COTTON CRAFT (PVT) LTD PAKISTAN

Plot #407, 408- Sundar Industrial Estate, Raiwind Road, Lahore
Ph : +92-42-35298604-06 , Email: info@cottoncraftpvt.com

Website : www.Cottoncraftpvt.com