

ABSORBENT RIBBON GAUZE (BPC)

DISCRIPTION:

Cotton cloth of plain weave bleached to a good white, clean, and reasonably free from weaving defects, cotton leaf and shells and having fast selvedge edges. It is colorless or almost odorless.it may be slightly off-white if it has been sterilized.

FEATURES:

- Sentry Medical Ribbon Gauzes are constructed from 100% cotton gauze, making them soft, gentle, and absorbent.
- The edges of this gauze is designed to eradicate loose fibers caused by fraying, making them ideal solutions for wound cavity packing
- Chemically inert and soft to give maximum protection and not cause irritation to skin.

INDICATIONS:

- Use to pack sinus, dental and throat cavities and to pack open infected wounds to assist healing.
- Used to stop the flow of blood absorb drainage.
- Use as temporary absorbent in open surgery and for packing large cavities produced by abdominal surgery.

PRECUATIONS:

- Sterile before use.
- Absorbent Ribbon Gauze swells in use and this should be taken into consideration when packing small cavities.

STORAGE:

- Keep in cool & dry place.
- Its absorbency may reduce considerably by medications, prolong storage.

AVAILABLE SIZES:

Gauze Ribbon 1.25cmx5M
Gauze Ribbon 1.25cmx90M
Gauze Ribbon 2.5cmx5M
Gauze Ribbon 2.5cmx90M



COTTON CRAFT (PVT) LTD PAKISTAN

Plot #407, 408- Sundar Industrial Estate, Raiwind Road, Lahore
Ph : +92-42-35298604-06 , Email: info@cottoncraftpvt.com

Website : www.Cottoncraftpvt.com