

OPTIC EYE PAD (BP)

DISCRIPTION:

It is an oval pad consisting of absorbent cotton, Faced on both sides with absorbent muslin cloth. It is sterile cotton pad individually packed in blister packing.

PRECUATIONS:

- Sterile before use.
- Before use after eye procedure ensures that the packet is intact.

FEATURES:

- These sterile Eye Pads are made with a soft non-woven cover, over absorbent cotton filler.
- Compressed edges help reduce linting.

STORAGE:

- Keep in cool & dry place.

INDICATIONS:

- Eye pads Protects eyes from dust and grime and protection from possible infection after surgical procedure.
- Eye pads protect injured eyes from further injury or prevent infection and strain.
- It can be safely used if there's a small cut to the eyelid.
- It is used to give comfort to the eye after dirt, smoke or girt exposure.

AVAILABLE SIZES:

6.0cmx8.0cm



COTTON CRAFT (PVT) LTD PAKISTAN

Plot #407, 408- Sundar Industrial Estate, Raiwind Road, Lahore
Ph.: +92-42-35298604-06 , Email: info@cottoncraftpvt.com

Website : Cottoncraftpvt.com